

# *PowerSoftMD* on Home Computer

It's OK with us if you install and run a back up copy of *PowerSoftMD* on your Home Computer. But, we do not support this directly. In other words, you are on your own setting it up and getting it to work. If you have problems please have a qualified computer technician assist you, not Data Tec.

1. Use the *PowerSoftMD* Update CD and pick the option from the CD Menu to "Add Work Station (new computer)". This will install the base *PowerSoftMD* files and insert the necessary support files in your Version of Windows.
2. Copy the entire **C:\EZW** folder and all subfolders from your main computer at the office and restore it in it's entirety to your home computer.

We recommend coping it to a USB device.

\*\*\* Be careful to use **Copy, NOT MOVE** \*\*\*

3. When you try to start *PowerSoftMD* (EasyMed) on the "Login" screen select the top toolbar "Remote Usage" option, then select "Run on Home Computer". This should make your home computer have a substituted drive that's the same letter as your office networked drive. Next, Login as normal.