

# *PowerSoftMD*

## Updates via the Internet

**WARNING** – We always highly recommend making sure you back up before applying any changes.

### **From a Work Station or the Main Computer.**

1. Do one of the following:
  - Open PowerSoftMD, from the Primary Menu select the Top toolbar “Tools” option then “Internet Download Updates”.
  - or
  - Using the Version Update Reminder Screen select the Green button.

You will be prompted with a message telling which **drive letter** you will Need to use when you unzip the download file, **make note of it.**

2. This should link to our update file and allow you to download it.
3. Once downloaded, change the folder name from **X:\** to **Z:\** where **Z** is the network drive letter you use to network. (Typically H:\ or G:\, whatever letter you made note of in step 1 above.)
4. Next click the Unzip button.
5. When prompted for a password it is **(Call Data Tec** for latest Password)
6. After the files are unzipped your main computer will be updated, you **must update each work station.** Use the top Toolbar Menu “Tools” option, then “Work Station Update” to update each work station.

NOTE: If the unzip password doesn't work, contact Data Tec and ask for the Current Internet Update Password.