

"New Work Station Installation" Steps

1. Set up the Network Drive Mapping on the Work Station

- A. Make sure you have the Main Computer location where PowerSoftMD resides mapped to the new computer exactly the same way all your other work stations are mapped. **All Work Stations** must use the **same Main Computer Drive Letter**, for example: H: or W:
- B. The drive must be mapped where you can see the **\EZW** folder in the **root** directory for the mapped drive. For example if you look at the **H:** drive you can see **H:\EZW**
- C. Make sure the Work Station has complete and total access rights to the mapped drive, i.e. Read, Write, Create, Delete, etc.

2. Install PowerSoftMD basic files on New Work station

- A. Download the file **NewComp.exe** from our web site at:
<http://www.powersoftmd.com/updates/NewComp.exe>
- B. Follow the on screen prompts as normal. Always, select "NEXT".

3. Start PowerSoftMD on the New Computer

- A. Using the American Flag "PowerSoftMD" Icon on the new work station desktop launch PowerSoftMD.
- B. You should be asked if you are installing a new work station reply **YES**
- C. A screen will be displayed, fill in the Network Main Computer Drive Letter and then click the "**Complete Installation**" button. This will open a window and you will see files being copied from your server to the new Work Station.
- D. Once step 3 is complete, launch PowerSoftMD again. On the Login screen select the top "**Tools**" option then "**Drive Selection**". Next, select the down arrow and pick the Network Main Computer Drive Letter. This will change PowerSoftMD's internal pointers to the proper drive letter.
- E. Launch PowerSoftMD one more time and you should be able to Login Normally.

