

PowerSoftMD Certified Enhancements 2011

12-31-2011

0016 You can set a new SOAP Note screen option "Prompt for Smoking Status Age >17". This will cause you to be reminded to gather Smoking Info for patients 18 years old or older, if you haven't already done so.

0015 If you are using the New Version of the Schedule, by default when you use the Check In feature on the schedule the Additional Demographics Window will be Displayed if the patient doesn't have all their Meaningful Use Demographics on their record. You can deactivate this option from the Utilities Options Menu.

12-03-2011

0014 The CPOE - Computerized Provider Order Entry system now available from the Top Window option on the Patient Profile & SOAP Note screens, in addition to the Primary Menu. Access from the Primary Menu lets you access orders for all patients, while launching it from a specific patient screen only accesses that patient's orders.

11-25-2011

0013 There is a new Patient Reminders Window that lets you enter and log your Patient Reminder contacts. It also lets you specify the methods the patient Wants to be contacted. You can access this Window from the MU Central Screen by selecting "View/Update Patient Reminders". The USPS Mail entry is automatically updated if you print Recall cards from the forms menu.

11-23-2011

0012 The Primary Menu Top EMR Icon defaults to opening the EMR Patient Profile Screen instead of the EMR Desktop. To open the EMR Desktop use the drop down Arrow to the right of the EMR Icon.

11-22-2011

0011 There is new option to Generate and test Hash Codes for files. Allowing you to determine if a file or files have been changed. From the EMR Desktop Screen select the top "Tools" option then "Hash Code - File Compare".

11-20-2011

0010 There is a new option to let you Encrypt and/or Decrypt Files. From the EMR Desktop screen select the top "Tools" option then "Encrypt/Decrypt Files".

0009 You can use a new MU Export screen to create and export Clinical Summary and/or complete CCR Records (Continuity Of Care Record. From the Primary Menu or from a Patient's Profile screen, select the top "Tools" option then "Export MU Files".

11-13-2011

0008 A new report is available to search and list the MU EMR Record Audit Log. From the Primary Menu select the drop down arrow to the right of the EMR Icon, then select "Reports". Next in the Meaningful Use box select the "Audit Log Report" option.

11-11-2011

0007 You can no longer Delete System User-ID's. But, you can In-Activate them.

11-10-2011

0006 Version 1 of the CPOE Computer Provider Order Entry software is available by

PowerSoftMD Certified Enhancements 2011

Using the "CPOE" button on the Primary Menu.

11-09-2011

0005 A new option allows you to update many patient demographics from the EMR Patient Profile screen. Use the top toolbar "Tools" option then "Meaningful Use Demographics".

0004 On the EMR Patient Profile screen you can run a BMI Calculator by selecting the top toolbar "Tools" then "BMI Calculator".

11-07-2011

0003 A new report is available to search EMR Records using Meaningful Use Criteria. From the Primary Menu use the Drop Down Arrow next to the top EMR Icon, and select Reports. Then select the report in the Green Meaningful Use window labeled "Patient Search by Criteria".

10-30-2011

0002 There is a new EMR Disease Registry window available for patients. You can enter and update many required MU Problems, Lab Results, Immunizations, etc. From the Patient Profile Screen select the "Disease Registry" tab, or from the Patient's Problem List Screen select the new Yellow "Disease Registry" button.

0001 There is a new EMR Smoking Status window available for patients. You can enter and update the required MU Smoking Status detail. From the Profile Screen click on the "Smoking Status" field to update it. You may also update it by using the "Soc/Fam History" button on the SOAP Note screen, then use the "Update" button on the Smoking Status area.