

Imaging Results

(Menu Objective #3)

OBJECTIVE:

Imaging results consisting of the image itself and any explanation or other accompanying information are accessible through CEHRT.

MEASURE:

More than 10 percent of all tests whose result is one or more images ordered by the EP during the EHR reporting period are accessible through CEHRT.

EXCLUSION:

Any EP who orders less than 100 tests whose result is an image during the EHR reporting period; or any EP who has no access to electronic imaging results at the start of the EHR reporting period.

Government Website for Imaging Results:


www.cms.gov/Regulations-and-Guidance/Legislation/EHRIncentivePrograms/downloads/Stage2_EPMenue_3_ImagingResults.pdf

Importing Image Results:

To Import Image Results, you need to record the imaging Orders, see instructions for CPOE.

The image files must be in a folder created for importing images from. Call us to set this up.


1) At the primary menu, click “Import Images” button in the bottom left corner.



Imaging Results

(Menu Objective #3)


2) Click on the image file you are importing



3) Windows will open the image selected, after verifying it's the right one, close it.




4) Click Import Screen, and click Yes to verify you closed the image.



Imaging Results

(Menu Objective #3)

5) Select the document folder to import to, then click Import now



6) Click the image order this result is for, and click Yes to mark it completed.

