

## **CLINICAL LAB-TEST RESULTS**

(Core Objective #10)

---

**OBJECTIVE:** Incorporate clinical lab-test results into Certified EHR Technology (CEHRT) as structured data.

**MEASURE:** More than 55 percent of all clinical lab tests results ordered by the EP during the EHR reporting period whose results are either in a positive/negative or numerical format are incorporated in Certified EHR Technology as structured data.

**EXCLUSION:** Any EP who orders no lab tests where results are either in a positive/negative affirmation or numeric format during the EHR reporting period.

---

**NUMERATOR:** (Lab-Tests Credited)

Number of matching Lab-Tests where the results are recorded as structured data.

**DENOMINATOR:** (Lab-Tests Considered)

Total number of Lab-Tests ordered by the provider whose results are numeric or positive/negative.

### **HOW TO:**

If you are receiving the Lab Results Electronically:

- 1) From the Primary Menu, click on the "Import HL7" button.
- 2) Select the "HL7" file Name.
- 3) View the Lab Result as an option.
- 4) Select the Date of Service and Import Folder.
- 5) Click the "Import HL7 Lab Results" button.

**Or** if you are scanning in or importing the Lab Results document (PDF).

- 1) From the Primary Menu, click on the "LabTrack: button.
- 2) On a lab request, click the date returned field and put the date of the return.
- 3) On the pop up, select import & import the HL7 results.

**Or** if you are using the Patient's Disease Registry.

- 1) From the Primary Menu use the EMR Drop down arrow and select "Disease Registry"
- 2) Select the Patient by Name.
- 3) On the Disease Registry, select the "Lab Results" tab.
- 4) Enter the results, then click the "Save/Exit" button.

1/20/2015