

MEDICATION RECONCILIATION

(Core Objective #14)

OBJECTIVE: The EP who receives a patient from another setting of care or provider of care or believes an encounter is relevant should perform medication reconciliation.

MEASURE: The EP who performs medication reconciliation for more than 50 percent of transitions of care in which the patient is transitioned into the care of the EP.

EXCLUSION: Any EP who was not the recipient of any transitions of care during the EHR reporting period.

NUMERATOR: (Patients Credited)

Number of matching patients who had medication reconciliation performed.

DENOMINATOR: (Patients Considered)

Total number of new patients transitioned to the provider in the Reporting Period.

HOW TO:

- 1) Import a CCDA from another EP to your computer system.
Use the eScript Interface "Secure Messaging" button to get the CCDA file.
- 2) From the patient's SOAP Note screen, select the top toolbar "MU Central" option.
- 3) Select the "Clinical Information Reconciliation" button.
- 4) Use the "Browse to Select New CCDA" button and click on the CCDA file name.
- 5) In the "List of Patient CCDA Files", click on the file name.
- 6) Use the "Allergies", "Medications", and "Problems" buttons to create a reconciled list for each category.
- 7) Exit the "CCDA Reconciliation" screen.
- 8) Open the patients eScripts screen, this will upload the reconciled Allergies and Medications.

Open the MU Central from the SOAP note, click on Clinical Information Reconciliation and reconcile the patients medications against a CCDA from their prior place of care.

1/20/2015